

Accessory

Disposable Plastic Dropper (In Bag)

Sample Diluent Buffer (In Bottle)

Cassette / Midstream (Empty Housing)

Sterile Safety Lancet (CE Mark)

Throat / Nasal Swab (CE Mark)

Sample Collection Tube (Secretion)

Alcohol Pad (CE Mark)

Aluminum Foil Pouch (Customized)


Biohazardous Waste Bag (Customized)

Sample Collection Tube (Stool)

Desiccant (CE Mark)

Pipette

Centrifuge



★ Remark : Accessories are optional and flexible for OEM requirements.
Equipments are optional for clients and / or end-users.

FAQs

? Basic Info (NewScen · Manufacturer)

📄 About us :

- MOQ : 200 Uncut Sheets
- Method : Lateral Flow Assay (Immunochromatography Assay)
- Capacity : $\geq 10,000$ Uncut Sheets (Regular Daily Output / Production Line)
- Category : 10 Categories (57 Kinds of Uncut Sheets in 10 Categories)

? Trial Order (OEM · Private Labeling)

📄 About you :

- Standard : QC Inspection Standard (Preferably in a Check Form, or the descriptions in details including requirements for the Performance / Appearance, etc.)
- Design : Design and Dimensions Drawing with precise sheet component position in details
- Format : Cassette / Midstream / Strip (Preferably in Picture with Dimension)


? Sample (Regular Format · Evaluation)

📄 For uncut sheet sample from NewScen, please contact us by website / e-mail / social media :

Website : www.newscentestkit.com

E-mail : export@newscen.com

Social Media : Facebook / Youtube / Linkedin / Instagram / Twitter



Contact NewScen Team :

Add : No. 65, 6th Street, Tianjin TEDA, 300457, Tianjin, P.R.China

Web : www.newscentestkit.com (Uncut Sheet)

Web : www.newscenbiotech.com (Rapid Test)

Email : export@newscen.com



Uncut Sheet




Rapid Test



Lateral Flow Assay UNCUT SHEET

2026 VERSION

PRODUCT CATALOG



Category of Uncut Sheet

INFECTIOUS DISEASE TEST KIT

	UNCUT SHEET	PRODUCT DESCRIPTION	SPECIMEN	SPECIFICATION
★	HIV 1+2	3 Lines : HIV (1+2) Antibody Rapid Test Kit	WHOLE BLOOD / SERUM / PLASMA	Cassette : 300 * 70 / 60 / 65 mm Strip / Dipstick : 300 * 80 / 70 mm
	HIV 1/2	2 Lines : HIV (1/2) Antibody Rapid Test Kit		
★	HCV	HCV Antibody Rapid Test Kit		
★	HBsAg	Hepatitis B Surface Antigen (HBsAg) Diagnostic Test Kit		
	HBsAb	Hepatitis B Surface Antibody (HBsAb) Diagnostic Test Kit		
	HEV	Hepatitis E Antibody (HEV IgM) Diagnostic Test Kit		
★	T P	Treponema Pallidum (TP) Antibody Rapid Test Kit		
	T B	M. Tuberculosis (TB IgG) Antibody Rapid Test Kit	SERUM / PLASMA (WHOLE BLOOD NOT APPLIED)	Cassette : 300 * 70 / 60 / 65 mm
	HAV	Hepatitis A Antibody (HAV IgM) Diagnostic Test Kit		
	HPAb	H. Pylori (HP) Antibody Rapid Test Kit	WHOLE BLOOD / SERUM / PLASMA	
★	HPAg	H. Pylori (HP) Antigen Rapid Test Kit	STOOL	
	ADENOVIRUS (Ag)	Adenovirus Antigen Rapid Test Kit		
	NORO (Ag)	Norovirus (NORO) Antigen Rapid Test Kit		
	ROTA (Ag)	Rotavirus (ROTA) Antigen Rapid Test Kit		
	TYPHOID (Ag)	Typhoid Antigen Rapid Test Kit		
	TYPHOID (IgG/IgM)	Typhoid IgG/IgM Antiobdy Rapid Test Kit	WHOLE BLOOD / SERUM	
	HSV (IgG)	Herpes Simplex Virus (HSV Type II) IgG Rapid Test Kit	SERUM / PLASMA (WHOLE BLOOD NOT APPLIED)	
	RSV (Ag)	Respiratory Syncytial Virus (RSV) Antigen Rapid Test Kit	NASOPHARYNGEAL SWAB	
	STREP A (Ag)	Group A Streptococcal Antigen (STREP A) Antigen Rapid Test Kit	SECRETION	
	CHLAMYDIA (Ag)	Chlamydia Antigen Rapid Test Kit		

CARDIAC MARKER TEST KIT

UNCUT SHEET	PRODUCT DESCRIPTION	SPECIMEN	SPECIFICATION
D-Dimer	D-Dimer Quantitative Rapid Test Kit	WHOLE BLOOD / PLASMA	Cassette : 300 * 70 / 60 / 65 mm
cTnl	Cardiac Troponin-I (cTnl) Quantitative Rapid Test Kit	WHOLE BLOOD / SERUM / PLASMA	

TROPICAL (MOSQUITO TRANSMITTED) DISEASE TEST KIT

	UNCUT SHEET	PRODUCT DESCRIPTION	SPECIMEN	SPECIFICATION
	MAR (Pf/Pan)	Malaria Pf/Pan Antigen Rapid Test Kit	WHOLE BLOOD (SERUM / PLASMA NOT APPLIED)	Cassette : 300 * 70 / 60 / 65 mm
★	DENGUE (IgG/IgM)	Dengue IgG/IgM Rapid Test Kit	WHOLE BLOOD / SERUM / PLASMA	
★	DENGUE (NS1 Ag)	Dengue NS1 Antigen Rapid Test Kit		
	ZIKA (IgG/IgM)	Zika Virus IgG/IgM Rapid Test Kit		
	ZIKA (NS1 Ag)	Zika NS1 Antigen Rapid Test Kit		
	CHIK (IgG/IgM)	Chikungunya IgG/IgM Rapid Test Kit		

DRUGS OF ABUSE (DOA) TEST KIT

UNCUT SHEET	PRODUCT DESCRIPTION	SPECIMEN	SPECIFICATION
5-IN-1	One Step (MOR / MET / KET / THC / MDMA) Multiple Drug Test Kit	URINE	Cassette : 300 * 70 / 60 / 65 mm Strip / Dipstick : 300 * 80 / 70 mm
★ MOR	Morphine (MOR) Gold Rapid Screen Test Kit		
★ MET	Methamphetamine (MET) Gold Rapid Screen Test Kit		
KET	Ketamine (KET) Rapid Test Kit		
THC	Tetrahydrocannabinol (THC) Rapid Test Kit		
MDMA	Methylenedioxymetham-Phetamine (MDMA) Rapid Test Kit		
AMP	Amphetamine (AMP) Rapid Test Kit		
BZD	Benzodiazepine (BZD) Rapid Test Kit		
COC	Cocaine (COC) Rapid Test Kit		
MTD	Methadone (MTD) Rapid Test Kit		

RESPIRATORY DISEASE (COVID-19 / INFLUENZA) TEST KIT

UNCUT SHEET	PRODUCT DESCRIPTION	SPECIMEN	SPECIFICATION
COVID-19 (Ag)	COVID-19 Antigen Rapid Test Cassette	THROAT / NASAL SWAB	Cassette : 300 * 70 / 60 / 65 mm
COVID-19 Ag & Flu (A+B)	COVID-19 Ag. & Influenza Combo Rapid Test Kit		
Flu (A/B)	Influenza (A/B) Rapid Test Kit		
Flu (A+B)	Influenza (A+B) Combo Rapid Test Kit		
COVID-19 (IgG/IgM)	3 Lines : COVID-19 IgG/IgM Rapid Test Cassette	WHOLE BLOOD / SERUM / PLASMA	
COVID-19 (NAb)	COVID-19 Neutralizing Antibody Rapid Test Cassette		

TUMOR MARKER TEST KIT

UNCUT SHEET	PRODUCT DESCRIPTION	SPECIMEN	SPECIFICATION
PSA	Prostate Specific Antigen (PSA) Diagnostic Test Kit	WHOLE BLOOD / SERUM / PLASMA	Cassette : 300 * 70 / 60 / 65 mm (FOB Strip : 300 * 80 / 70 mm)
AFP	Alpha Fetoprotein (AFP) Quantitative Rapid Test Kit		
CEA	Carcino-Embryonic Antigen Quantitative Rapid Test Kit		
★ FOB	One-Step Fecal Occult Blood (FOB) Diagnostic Test Kit	STOOL	
T F	Transferrin (TF) Rapid Test Kit		

HORMONE TEST KIT

UNCUT SHEET	PRODUCT DESCRIPTION	SPECIMEN	SPECIFICATION
HCG	Human Chorionic Gonadotropin (hCG) Tests (Pregnancy tests)	URINE	Midstream : 300 * 70 / 60 / 65 mm
L H	Luteinizing Hormone (LH) Tests (Ovulation tests)		Cassette : 300 * 70 / 60 / 65 mm
FSH	Follicle Stimulating Hormone (FSH) Tests (Menopause Tests)		
TSH	Human Thyroid Stimulating Hormone (TSH) Quantitative Rapid Test Kit	WHOLE BLOOD / SERUM / PLASMA	Cassette : 300 * 70 / 60 / 65 mm

IRON DEFICIENCY TEST KIT

UNCUT SHEET	PRODUCT DESCRIPTION	SPECIMEN	SPECIFICATION
★ S F	Serum Ferritin (SF) Rapid Test Kit	WHOLE BLOOD / SERUM / PLASMA	Cassette : 300 * 70 / 60 / 65 mm

INFLAMMATORY MARKER TEST KIT

UNCUT SHEET	PRODUCT DESCRIPTION	SPECIMEN	SPECIFICATION
CRP	C-Reactive Protein (CRP) Quantitative Rapid Test Kit	WHOLE BLOOD / SERUM / PLASMA	Cassette : 300 * 70 / 60 / 65 mm
PCT	Human Procalcitonin (PCT) Quantitative Rapid Test Kit		

RENAL INJURY TEST KIT

UNCUT SHEET	PRODUCT DESCRIPTION	SPECIMEN	SPECIFICATION
mAlb	Urine Microalbumin (mAlb) Quantitative Rapid Test Kit	URINE	Cassette : 300 * 70 / 60 / 65 mm

★ Star Products : Hot Sale Uncut Sheet